

## Product Datasheet

# Human Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit (orb1665369)

<b>Catalog Number</b>	orb1665369
<b>Category</b>	Assays and Kits
<b>Description</b>	Human Interleukin 8 Receptor Beta (IL8Rb) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detect Range: 0.781-50 ng/mL. Sensitivity: 0.25 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)