

Product Datasheet

Rat Interleukin 18 Binding Protein (IL18BP) ELISA Kit (orb1665537)

Catalog Number	orb1665537
Category	Assays and Kits
Description	Rat Interleukin 18 Binding Protein (IL18BP) ELISA Kit
Reactivity	Rat
Tested applications	ELISA
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.215 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)