

Product Datasheet

Mouse Insulin Like Growth Factor 1 (IGF1) ELISA Kit (orb1665676)

Catalog Number	orb1665676
Category	Assays and Kits
Description	Mouse Insulin Like Growth Factor 1 (IGF1) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 29 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)