

Product Datasheet

Human Isocitrate Dehydrogenase 2, mitochondrial (IDH2) ELISA Kit (orb1665723)

Catalog Number	orb1665723
Category	Assays and Kits
Description	Human Isocitrate Dehydrogenase 2, mitochondrial (IDH2) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)