

Product Datasheet

Mouse Growth Regulated Oncogene Beta (GROb) ELISA Kit (orb1666001)

Catalog Number	orb1666001
Category	Assays and Kits
Description	Mouse Growth Regulated Oncogene Beta (GROb) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.0 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)