

Product Datasheet

Rat Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit (orb1666074)

Catalog Number	orb1666074
Category	Assays and Kits
Description	Rat Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit
Reactivity	Rat
Tested applications	ELISA
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.05 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road

Cambridge

CB5 8LA

United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive

Research Triangle Park

Durham

NC 27713

United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)