

Product Datasheet

Rat Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit (orb1666074)

Catalog Number	orb1666074
Category	Assays and Kits
Description	Rat Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit
Reactivity	Rat
Tested applications	ELISA
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.05 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)