

Product Datasheet

Mouse Gastric Inhibitory Polypeptide (GIP) ELISA Kit (orb1666130)

Catalog Number	orb1666130
Category	Assays and Kits
Description	Mouse Gastric Inhibitory Polypeptide (GIP) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 24.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)