

Product Datasheet

Mouse Estrogen Receptor Alpha (ERα) ELISA Kit (orb1666574)

Catalog Number	orb1666574
Category	Assays and Kits
Description	Mouse Estrogen Receptor Alpha (ERα) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.146 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)