

Product Datasheet

Mouse Epidermal Growth Factor Receptor 2 (EGFR2) ELISA Kit (orb1666646)

Catalog Number	orb1666646
Category	Assays and Kits
Description	Mouse Epidermal Growth Factor Receptor 2 (EGFR2) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.207 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)