

## Product Datasheet

# Human Cytochrome b-245 Beta Polypeptide (CYBb) ELISA Kit (orb1666884)

|                            |                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1666884                                                                                                                                                                                                  |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                                             |
| <b>Description</b>         | Human Cytochrome b-245 Beta Polypeptide (CYBb) ELISA Kit                                                                                                                                                    |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                       |
| <b>Tested applications</b> | ELISA                                                                                                                                                                                                       |
| <b>Application notes</b>   | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.055 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: tissue homogenates, cell lysates or other biological fluids. |
| <b>Note</b>                | For research use only                                                                                                                                                                                       |

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)