

Product Datasheet

Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit (orb1667335)

Catalog Number	orb1667335
Category	Assays and Kits
Description	Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit
Reactivity	Mouse
Tested applications	ELISA
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.032 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)