

Product Datasheet

Human Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit (orb1667336)

Catalog Number	orb1667336
Category	Assays and Kits
Description	Human Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 33 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)