

Product Datasheet

Human Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit (orb1667974)

| | |
|----------------------------|---|
| Catalog Number | orb1667974 |
| Category | Assays and Kits |
| Description | Human Alkaline Sphingomyelinase (Alk-Smase) ELISA Kit |
| Reactivity | Human |
| Tested applications | ELISA |
| Application notes | Detect Range: 3.12-200 ng/mL. Sensitivity: 1.39 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates or other biological fluids. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)