

Product Datasheet

Sheep Ghrelin/obestatin prepropeptide (GHRL) ELISA Kit (orb1668573)

Catalog Number	orb1668573
Category	Assays and Kits
Description	Sheep Ghrelin/obestatin prepropeptide (GHRL) ELISA Kit
Reactivity	Sheep
Tested applications	ELISA
Application notes	Detect Range: 123.46-10000 pg/mL. Sensitivity: 48.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)