

Product Datasheet

Rabbit Cross linked C-telopeptide of type II collagen (CTX-II) ELISA Kit (orb1669598)

Catalog Number	orb1669598
Category	Assays and Kits
Description	Rabbit Cross linked C-telopeptide of type II collagen (CTX-II) ELISA Kit
Reactivity	Rabbit
Tested applications	ELISA
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 40 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)