

Product Datasheet

Guinea pig Interleukin-1 beta (IL1B) ELISA Kit (orb1670425)

Catalog Number	orb1670425
Category	Assays and Kits
Description	Guinea pig Interleukin-1 beta (IL1B) ELISA Kit
Reactivity	Guinea pig
Tested applications	ELISA
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 15.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)