

Product Datasheet

Human Lecithin Cholesterol Acyltransferase (LCAT) ELISA Kit (orb1670582)

Catalog Number	orb1670582
Category	Assays and Kits
Description	Human Lecithin Cholesterol Acyltransferase (LCAT) ELISA Kit
Reactivity	Human
Tested applications	ELISA
Application notes	Detect Range: 0.469-30 ng/mL. Sensitivity: 0.19 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)