

Product Datasheet

Sheep Luteinizing Hormone(LH)ELISA Kit (orb1670825)

Catalog Number	orb1670825
Category	Assays and Kits
Description	Sheep Luteinizing Hormone(LH)ELISA Kit
Reactivity	Sheep
Tested applications	ELISA
Application notes	Detect Range: 2 mIU/mL-75 mIU/mL. Sensitivity: 0.5 mIU/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)