

Product Datasheet

Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit (orb1670869)

Catalog Number	orb1670869
Category	Assays and Kits
Description	Pig Tumor necrosis factor (TNF/TNFA/TNFSF2) ELISA kit
Reactivity	Porcine
Tested applications	ELISA
Application notes	Detect Range: 0.047 ng/mL-3 ng/mL. Sensitivity: 0.040 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)