

## Product Datasheet

# Tumor Necrosis Factor (TNF) alpha Recombinant Protein (orb1671197)


|                            |                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1671197                                                                                                                                          |
| <b>Description</b>         | Tumor Necrosis Factor (TNF) alpha Recombinant Protein                                                                                               |
| <b>Tested Applications</b> | CLIA, ELISA, LF, WB                                                                                                                                 |
| <b>Preservatives</b>       | PBS, pH 8.0, and 8 M Urea containing 0.1% ProClin300.                                                                                               |
| <b>Form/Appearance</b>     | liquid                                                                                                                                              |
| <b>Concentration</b>       | batch dependent                                                                                                                                     |
| <b>Storage</b>             | Upon receipt, store at 2-8°C                                                                                                                        |
| <b>Note</b>                | For research use only                                                                                                                               |
| <b>Application notes</b>   | Suitable for use in ELISA, CLIA, Lateral Flow, Western Blot, and other immune tests                                                                 |
| <b>Purity</b>              | > 90% pure by SDS-PAGE. Ligand affinity purified.                                                                                                   |
| <b>Source</b>              | E. coli                                                                                                                                             |
| <b>MW</b>                  | This recombinant protein contains full length of Human TNF-alpha with a MW of 18.6 kDa. It is recognized by an anti-human TNF-alpha in ELISA tests. |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                      |

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



SDS-Page Image

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)