

Product Datasheet

Beta 2 adrenergic receptor Antibody [R11/E1], Mouse IgG1 (orb1671588)

Catalog Number	orb1671588
Description	Beta 2 adrenergic receptor Antibody [R11/E1], Mouse IgG1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Immunogen	This antibody was raised by immunizing mice with beta-2 adrenergic receptor.
Target	ADRB2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Clone Number	R11/E1
Antibody Type	Primary Antibody
Uniprot ID	P07550
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)