

Product Datasheet

CD19 Antibody [FMC63], Mouse IgG2a (orb1671950)


Description	CD19 Antibody [FMC63], Mouse IgG2a
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human polymphocytic leukaemia cell line JVM3.
Target	CD19
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Clone Number	FMC63
Antibody Type	Primary Antibody
Uniprot ID	P15391
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Immunohistochemical staining of rat spleen tissue using anti-CD19 antibody (FMC63). Anti-CD19 staining of formaldehyde fixed paraffin embedded rat spleen tissue - at 40x magnification. The mouse IgG 2A-chimeric version of FMC63 (orb1671950) was used to stain samples at a concentration of 3 ug/ml.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)