

Product Datasheet

CD16/CD32 Antibody [2.4G2], Rabbit IgG (orb1672261)

Description	CD16/CD32 Antibody [2.4G2], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IF, IHC, IP
Immunogen	Mouse BALB/c Macrophage J774 Cell Line.
Target	Fcgr2b
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Clone Number	2.4G2
Antibody Type	Primary Antibody
Uniprot ID	P08508, P08101
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)