

Product Datasheet

MYC Antibody (orb1672331)

Catalog Number	orb1672331
Category	Antibodies
Description	MYC Antibody
Target	MYC
Clonality	Recombinant
Isotype	IgG1 kappa
Conjugation	Unconjugated
Reactivity	Human
Concentration	batch dependent
Buffer/Preservatives	PBS with 0.02% Proclin 300.
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
UniProt ID	P01106
Tested applications	IF, IHC-Fr, IHC-P, IP, WB
Antibody Type	Primary Antibody
Clone Number	9E10
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)