

Product Datasheet

c-myc epitope tag Antibody [9E10], Mouse IgG2c (orb1672332)

Description	c-myc epitope tag Antibody [9E10], Mouse IgG2c
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Target	MYC
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2c kappa
Clonality	Recombinant
Clone Number	9E10
Antibody Type	Primary Antibody
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)