

Product Datasheet

Agarose Immobilized IgG Fraction of anti-Goat IgG [H&L] [Rabbit] (orb1674980)

Description	Agarose Immobilized IgG Fraction of anti-Goat IgG [H&L] [Rabbit]
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Unconjugated
Tested Applications	Purification
Immunogen	Goat IgG whole molecule
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Storage	Store vial at 4° C prior to opening. DO NOT FREEZE. Expiration date is one (1) year from date of opening.
Note	For research use only
Application notes	Not specified. Optimal titers for applications should be determined by the researcher.
Antibody Type	Secondary Antibody
Purity	This product is an IgG fraction coupled to activated agarose. Sufficient capacity is provided to bind a minimum of 5 mg of pure Goat IgG.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)