

Product Datasheet

Fluorescein Conjugated IgG fraction of anti-Rabbit IgG [H&L] Mouse] (orb1674985)

Catalog Number	orb1674985
Description	Fluorescein Conjugated IgG fraction of anti-Rabbit IgG [H&L] Mouse]
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	FACS
Immunogen	Rabbit IgG whole molecule
Preservatives	0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4
Concentration	10.0 mg/ml (by UV absorbance at 280 nm)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Not specified. Optimal titers for applications should be determined by the researcher
Antibody Type	Secondary Antibody
Purity	This product is an IgG fraction purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Mouse Serum, Rabbit IgG and Rabbit Serum.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)