

Product Datasheet

MAP2 monoclonal antibody (orb1677294)

Description	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% sodium azide.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human MAP2. The exact sequence is proprietary.
Preservatives	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% sodium azide.
Form/Appearance	Mouse IgG1 kappa. Supplied in crude ascites with 0.01% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Antibody Type	Primary Antibody
MW	~ 59 kDa
Uniprot ID	P11137
Dilution Range	WB (1/500 - 1/1000), IHC (1/50 - 1/200), IF/ICC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)