

Product Datasheet

LXR agonist 2 (orb1685892)

Catalog Number	orb1685892
Category	Small Molecules
Description	LXR agonist 2 is a potent small molecule activator of Liver X Receptor (LXR) that functions by stabilizing the NCOA1 coactivator. This tool compound is used in research to study LXR-mediated pathways in lipid metabolism and inflammation, applicable in both cellular and animal models.
Target	Liver X Receptor,Others
MW	586.16
Formula	$C_{35}H_{40}ClN_3O_3$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)