


Product Datasheet

Toreforant (orb1706521)

| | |
|------------------------|--------------------------------------------------------------------------------------------------------------------|
| Catalog Number | orb1706521 |
| Category | Small Molecules |
| Description | Toreforant is an effective and selective histamine H4 receptor antagonist (Ki at the human receptor: 8.4 nM). |
| Target | Histamine Receptor |
| MW | 392.54 |
| CAS Number | 952494-46-1 |
| Formula | C ₂₃ H ₃₂ N ₆ |
| SMILES | CN1CCC(CCCNc2ncc(-c3nc4cc(C)cc(C)c4[nH]3)c(C)n2)CC1 |
| Storage | Powder: -20°C for 3 years In solvent: -80°C for 1 year Shipping with blue ice/Shipping at ambient temperature. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |



Chemical structure of toreforant

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)