


Product Datasheet

Toreforant (orb1706521)

Catalog Number	orb1706521
Category	Small Molecules
Description	Toreforant is an effective and selective histamine H4 receptor antagonist (Ki at the human receptor: 8.4 nM).
Target	Histamine Receptor
MW	392.54
CAS Number	952494-46-1
Formula	C ₂₃ H ₃₂ N ₆
SMILES	CN1CCC(CCCNc2ncc(-c3nc4cc(C)cc(C)c4[nH]3)c(C)n2)CC1
Storage	Powder: -20°C for 3 years In solvent: -80°C for 1 year Shipping with blue ice/Shipping at ambient temperature.
Note	For research use only
Expiration Date	12 months from date of receipt.



Chemical structure of toreforant

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)