

Product Datasheet

S1P1 Agonist III (orb1706670)

Catalog Number	orb1706670
Category	Small Molecules
Description	S1P1 Agonist III
Target	LPL Receptor,Others
Purity	99.08%
MW	415.37
Target Areas	S1P Receptor
Solubility (25°C)	10% DMSO+40% PEG300+5% Tween-80+45% Saline:2.5 mg/mL (6.02 mM);DMSO:48 mg/mL (115.56 mM)
CAS Number	1324003-64-6
Formula	C ₂₁ H ₁₆ F ₃ N ₃ O ₃
SMILES	O=C(NC(=O)C=1C=NC=CC1OC)NC=2C=CC(C=3C=CC=CC3)=C(C2)C(F)(F)F
Storage	store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year Shipping with blue ice/Shipping at ambient temperature.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)