

Product Datasheet

Ac4GalNAz (orb1710785)

Catalog Number	orb1710785
Category	Small Molecules
Description	Ac4GalNAz is an alkyl chain-based PROTAC linker utilized in PROTAC synthesis.
Target	PROTAC Linker
Purity	99.95%
MW	430.37
Target Areas	PROTAC Linker
Solubility (25°C)	DMSO:10 mM;10% DMSO+40% PEG300+5% Tween 80+45% Saline:2 mg/mL (4.65 mM)
CAS Number	653600-56-7
Formula	C ₁₆ H ₂₂ N ₄ O ₁₀
SMILES	<chem>CC(=O)OC[C@H]1OC(OC(C)=O)[C@H](NC(=O)CN=[N+]=[N-])[C@@H](OC(C)=O)[C@H]1OC(C)=O</chem>
Storage	Powder: -20°C for 3 years In solvent: -80°C for 1 year Shipping with blue ice/Shipping at ambient temperature.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558