

## Product Datasheet

### Rabbit anti-Thrombomodulin Recombinant Monoclonal Antibody (orb1710820)

|                             |                                                                                                                      |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb1710820                                                                                                           |
| <b>Category</b>             | Antibodies                                                                                                           |
| <b>Description</b>          | Thrombomodulin Antibody                                                                                              |
| <b>Target</b>               | Thrombomodulin                                                                                                       |
| <b>Clonality</b>            | Recombinant                                                                                                          |
| <b>Species/Host</b>         | Rabbit                                                                                                               |
| <b>Conjugation</b>          | Unconjugated                                                                                                         |
| <b>Reactivity</b>           | Mouse                                                                                                                |
| <b>Form/Appearance</b>      | Whole IgG                                                                                                            |
| <b>Concentration</b>        | 1000 µg/ml                                                                                                           |
| <b>Buffer/Preservatives</b> | Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free                                           |
| <b>Immunogen</b>            | Residues 17-517 (ECD)                                                                                                |
| <b>UniProt ID</b>           | <b>P15306</b>                                                                                                        |
| <b>Tested applications</b>  | ICC, IHC, IP, WB                                                                                                     |
| <b>Dilution range</b>       | WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF: |
| <b>Application notes</b>    | Format: Whole IgG                                                                                                    |
| <b>Antibody Type</b>        | Primary Antibody                                                                                                     |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

|                        |                                                                                                                                         |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clone Number</b>    | BLR260L                                                                                                                                 |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>            | For research use only                                                                                                                   |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                         |

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)