

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) Antibody (orb1711813)

Catalog Number	orb1711813
Category	Antibodies
Description	TAG-72 / CA72.4 (Tumor-Associated Glycoprotein); Clone SPM536 (Concentrate)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1, kappa.
Conjugation	Unconjugated
Reactivity	Bovine, Canine, Human, Rat
Immunogen	Purified TAG-72 protein.
Application notes	Jurkat cells. Breast or lung carcinoma.
Specificity	Recognizes an oncofetal antigen of 220kDa, identified as a tumor-associated glycoprotein (TAG-72) with properties of a mucin.
Clone Number	SPM536.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)