

Product Datasheet

MALT1 (MALT-Lymphoma Marker) Antibody (orb1712270)

Catalog Number	orb1712270
Category	Antibodies
Description	MALT1 (MALT-Lymphoma Marker); Clone MT1/410 (Concentrate)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1, kappa
Conjugation	Unconjugated
Reactivity	Human
Immunogen	Human MALT1 recombinant fragment (aa 701-808)
Application notes	Jurkat, Daudi, or HeLa cells. Tonsil or lymphoma.
Specificity	Highly expressed in peripheral blood mononuclear cells. Detected at lower levels in bone marrow, thymus and lymph node, and at very low levels in colon and lung.
Clone Number	MT1/410
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)