

Product Datasheet

CD63 (Late Endosomes Marker) Antibody (orb1712528)

Catalog Number	orb1712528
Description	CD63 (Late Endosomes Marker); Clone MX-49.129.5 (Concentrate)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Immunogen	Full length CD63 of human origin
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	SK-MEL-28, HL60, THP-1 or NIH/3T3 cells. Melanoma or lymphoma.
Isotype	IgG1, kappa
Clonality	Monoclonal
Clone Number	MX-49.129.5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)