

Product Datasheet

Recombinant Human USP10 (orb1714836)

Catalog Number	orb1714836
Category	Proteins
Description	Recombinant human USP10 protein is expressed in E. coli system with purity greater than 80% as determined by SDS-PAGE.
Tag	N-terminal His-IF2DI Tag
Reactivity	Human
Buffer/Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM HEPES, 500 mM NaCl with 5% trehalose, pH 7.4
Purity	Greater than 80% as determined by SDS-PAGE
UniProt ID	Q14694
MW	44.1 kDa
Tested applications	ELISA, WB
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Source	E.Coli
Expression Region	1-219
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)