

Product Datasheet

14-3-3 sigma / Stratifin Immunizing Peptide (orb17237)

Catalog Number	orb17237
Description	SFN Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	This peptide is used as immunogen to produce SFN antibody (orb19148)
Protein Sequence	DLHTLSEDSYKDST
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)