

Product Datasheet

Plasmid DNA Extraction and Purification Antibody (orb1726449)

Catalog Number	orb1726449
Category	Assays and Kits
Description	Plasmid DNA Extraction and Purification Antibody
Storage	Mixed storage temperatures. Magnetic beads, and RNase A and Neutralization solutions stored at 2-8°C. All other reagents stored at room temperature.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)