

Product Datasheet

Agarose Gel DNA Extraction and Purification Antibody (orb1726453)

Catalog Number	orb1726453
Category	Assays and Kits
Description	Agarose Gel DNA Extraction and Purification Antibody
Storage	Mixed storage temperatures. Magnetic beads stored at 2-8°C. Other reagents stored at room temperature.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)