

Product Datasheet

Human Galectin-8 ELISA Kit (orb1728073)

Catalog Number	orb1728073
Category	Assays and Kits
Description	Human Galectin-8 ELISA Kit
Target	Galectin-8
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-W317
UniProt ID	O00214
Tested applications	ELISA
Specificity	Natural and recombinant human LGALS8
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	46.9 pg/ml - 3,000 pg/ml
Sensitivity	25 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Human Galectin-8 ELISA Kit

Human Galectin-8 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)