

Product Datasheet

Goat anti Rabbit IgG-HRP conjugated (orb1729599)

Description	Goat polyclonal secondary antibody raised against Rabbit IgG-Fc heavy chain and conjugated with HRP
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Rabbit IgG-Fc-Heavy chain
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Clone Number	NA
MW	NA
Hazard Information	For Research use only
Dilution Range	WB- 1:3000 - 1:5000 ELISA- 1:10000-1:20000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)