

Product Datasheet

Total RNA Purification Kit - Column Kit (orb1733629)

Catalog Number	orb1733629
Category	Assays and Kits
Description	Total RNA Purification Kit is designed for rapid, high purity and high yield isolation of total RNA from small amounts of various samples including blood, animal and plant tissue, bacteria and viruses. The spin column based method allows complete removal of inhibitors such as divalent cations and proteins. Due to elimination of phenol, handling of the kit is safe and no harmful waste is produced. The purified total RNA can be used in a number of downstream applications. The kit allows the purification of up to 100 µg RNA per preparation.
Storage	store at ambient temperature
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)