

Product Datasheet

Yeast DNA Preparation - Solution Kit (orb1733630)

| | |
|------------------------|---|
| Catalog Number | orb1733630 |
| Category | Assays and Kits |
| Description | Yeast DNA Preparation Kit is designed for convenient and fast isolation of genomic DNA from yeast cells. The solution based system minimizes DNA fragmentation that may be problematic in other spin-column/filtration based method. Because phenol or chloroform is not used it is safe and does not produce any harmful waste. Solution based genomic DNA purification kits guarantee minimal DNA fragmentation and yield DNA sized up to 150 kb. |
| Storage | store at ambient temperature |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)