

Product Datasheet

Random Octamers (orb1733641)

Catalog Number	orb1733641
Category	Tools
Description	Random Octamers are a mixture of oligonucleotides representing all possible sequences for an octamer. Random Octamers are used in DNA labeling by PCR (DOP-PCR) or cDNA synthesis by RT-PCR.
Form/Appearance	clear aqueous solution
Concentration	100 µM
Purity	by HPLC
MW	Theoretical MW: 2410 g/mol
Storage	store at -20 °C. avoid freeze/thaw cycles
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558