

Product Datasheet

SNRPA / U1A Antibody (clone 3F9-1F7) (orb1735051)

Catalog Number	orb1735051
Category	Antibodies
Description	SNRPA / U1A Antibody (clone 3F9-1F7)
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1,k
Conjugation	Unconjugated
Reactivity	Human
Concentration	0.5 mg/ml
Buffer/Preservatives	PBS, pH 7.4
Purification	Purified from ascites by Protein A
Immunogen	SNRPA (AAH00405, 1 aa ~ 283 aa) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Tested applications	ELISA, IF, IHC, IHC-P, WB
Dilution range	ELISA, IF (10 µg/ml), IHC, IHC-P (5 µg/ml), WB
Application notes	Further information: Western Blot (Cell lysate) - positive control HeLa. Western Blot (Transfected lysate) - positive control transfected 293T cell line. Immunofluorescence (10 ug/ml) - positive control HeLa cells. Sandwich ELISA (Recombinant protein).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	3F9-1F7
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.



Immunofluorescence of monoclonal antibody to SNRPA on HeLa cell (antibody concentration 10 ug/ml).



Western blot of SNRPA expression in transfected 293T cell line by SNRPA monoclonal antibody clone 3F9-1F7.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Detection limit for recombinant GST tagged SNRPA is approximately 0.03 ng/ml as a capture antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)