

Product Datasheet

Rat IgG1 kappa (Isotype control) Antibody (orb1736005)

Catalog Number	orb1736005
Category	Antibodies
Description	Rat IgG1 kappa (Isotype control) Antibody
Species/Host	Rat
Isotype	IgG1 kappa
Conjugation	Unconjugated
Concentration	1 mg/mL
Buffer/Preservatives	PBS/0,1%NaN3
Purification	Immunoaffinity chromatography (Purification on Mouse anti-Rat column)
Tested applications	ELISA, IA
Clone Number	IR871
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)