

Product Datasheet

PCR Master Mix(Red) (orb1736815)

Catalog Number orb1736815

Category Tools

Description . . PCR Master Mix is a ready-to-use conventional PCR master mix solution containing Taq DNA Polymerase, dNTP mix, MgCl₂, and optimized buffer system. The amplification can be performed by adding primers and templates during the reaction, which greatly simplifies the operation steps of the experiment. The product contains cresol red, and the PCR product can be directly electrophoresed. This product contains excellent stabilizer and can be stored for 3 months at 4°C. PCR products have 3'-dA overhangs for easy cloning into T-vectors

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)